



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Protein A Magnetic Beads IP/Co-IP Kit ROC-KBA-PA00-18

Article Name	Protein A Magnetic Beads IP/Co-IP Kit
Biozol Catalog Number	ROC-KBA-PA00-18
Supplier Catalog Number	KBA-PA00-18
Alternative Catalog Number	ROC-KBA-PA00-18
Manufacturer	Rockland Immunochemicals
Category	Kits/Assays
Application	IP
Concentration	5.0 mg/mL
Form	Suspension of Magnetic Beads
Application Dilute	IP: User Optimized, WB: User Optimized
Application Notes	TrueBlot Protein A magnetic beads IP/Co-IP kit can be used for immunoprecipitation (IP) and Co-IP experiments as well as for antibody purification. For IP, target specific antibody is incubated with cell lysate. Protein A magnetic beads are then incubat