



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Protein G Magnetic Beads IP/Co-IP Kit ROC-KBA-PG00-18

Article Name	Protein G Magnetic Beads IP/Co-IP Kit
Biozol Catalog Number	ROC-KBA-PG00-18
Supplier Catalog Number	KBA-PG00-18
Alternative Catalog Number	ROC-KBA-PG00-18
Manufacturer	Rockland Immunochemicals
Category	Kits/Assays
Application	IP
Concentration	5.0mg/mL by UV absorbance at 280 nm
Form	Suspension of Magnetic Beads
Application Dilute	IP: User Optimized, WB: User Optimized
Application Notes	TrueBlot Protein G magnetic beads IP/Co-IP kit can be used for immunoprecipitation (IP) and Co-IP experiments as well as for antibody purification. For IP, target specific antibody is incubated with cell lysate. Protein G magnetic beads are then incubat