



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## **Transcription Factor Targeted ChIP-seq Kit ROC-KOA0882**

Article Name	Transcription Factor Targeted ChIP-seq Kit
Biozol Catalog Number	ROC-KOA0882
Supplier Catalog Number	KOA0882
Alternative Catalog Number	ROC-KOA0882
Manufacturer	Rockland Immunochemicals
Category	Kits/Assays
Application Dilute	ChIP: User Optimized
Application Notes	The different steps of the ChIP assay are cell fixation (cross-linking), chromatin shearing, immunoprecipitation, reverse cross-linking followed by DNA purification and analysis of the immunoprecipitated DNA. Suitable for 10 reactions.