



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

2X SDS-PAGE Sample Buffer without DTT or b-ME ROC-MB-018

Article Name	2X SDS-PAGE Sample Buffer without DTT or b-ME
Biozol Catalog Number	ROC-MB-018
Supplier Catalog Number	MB-018
Alternative Catalog Number	ROC-MB-018
Manufacturer	Rockland Immunochemicals
Category	Sonstiges
Application	SDS-PAGE
Conjugation	Unconjugated
Concentration	2X
Buffer	See application note.
Form	Liquid (sterile filtered)
Application Notes	2X SDS-PAGE Sample Loading Buffer is suitable for laboratory assays involved in protein biochemistry. 2X SDS-PAGE Sample Buffer is a concentrated stock solution and should be diluted appropriately with distilled, deionized water or equivalent to its fina