

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

ELISA Microwell Blocking Buffer with Stabilizer (Azide and Mercury Free) ROC-MB-064-1000

Article Name	ELISA Microwell Blocking Buffer with Stabilizer (Azide and Mercury Free)
Biozol Catalog Number	ROC-MB-064-1000
Supplier Catalog Number	MB-064-1000
Alternative Catalog Number	ROC-MB-064-1000
Manufacturer	Rockland Immunochemicals
Category	Kits/Assays
Application	ELISA
Conjugation	Unconjugated
Concentration	1X
Buffer	See application note.
Form	Liquid (sterile filtered)
Application Dilute	ELISA: User Optimized
Application Notes	This product is a ready-to-use 1X solution for stabilizing ELISA plates for storage prior to ELISA. After coating with antigen or antibody, wash the contents of the ELISA microwells and add a sufficient volume of ELISA Microwell Blocking Buffer with Stab