

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Diamond DNA Ladder 50-150bp ROC-MB-206-0500

| | |
|----------------------------|--|
| Article Name | Diamond DNA Ladder 50-150bp |
| Biozol Catalog Number | ROC-MB-206-0500 |
| Supplier Catalog Number | MB-206-0500 |
| Alternative Catalog Number | ROC-MB-206-0500 |
| Manufacturer | Rockland Immunochemicals |
| Category | Molekularbiologie |
| Application | SDS-PAGE |
| Conjugation | Unconjugated |
| Concentration | 112 µg/mL Sufficient to run approximately 100 tests |
| Buffer | See application note. |
| Form | Liquid (in 1x Loading Buffer) |
| Application Notes | Diamond Ready-to-Use DNA Ladder 50-150bp has a suggested load volume is 5 µL/well. |