



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

TrueBlot Streptavidin Magnetic Beads, SA ROC-S000-18-2

Article Name	TrueBlot Streptavidin Magnetic Beads, SA
Biozol Catalog Number	ROC-S000-18-2
Supplier Catalog Number	S000-18-2
Alternative Catalog Number	ROC-S000-18-2
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Conjugation	SA
Concentration	5.0 mg/mL
Buffer	0.1M Tris buffered saline, pH8.2
Form	Suspension of Magnetic Beads
Application Dilute	IP: User Optimized, WB: User Optimized
Application Notes	TrueBlot Streptavidin magnetic beads can be used for immunoprecipitation (IP) of antigens (using biotinylated antibodies) from a wide variety of sources and to Co-immunoprecipitate interaction complexes using biotinylated antibodies. Streptavidin magnet