

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## **TrueBlot Streptavidin Magnetic Beads, SA ROC-S000-18-5**

Article Name	TrueBlot Streptavidin Magnetic Beads, SA
Biozol Catalog Number	ROC-S000-18-5
Supplier Catalog Number	S000-18-5
Alternative Catalog Number	ROC-S000-18-5
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Conjugation	SA
Concentration	5.0 mg/mL
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Form	Suspension of Magnetic Beads
Application Dilute	IP: User Optimized, WB: User Optimized
Application Notes	TrueBlot Streptavidin magnetic beads are stored in phosphate buffered saline with detergent. TrueBlot Streptavidin magnetic beads can be used for immunoprecipitation (IP) of antigens (using biotinylated antibodies) from a wide variety of sources and to