



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

LNCaP Cell Nuclear Extract ROC-W09-001-GK1

Article Name	LNCaP Cell Nuclear Extract
Biozol Catalog Number	ROC-W09-001-GK1
Supplier Catalog Number	W09-001-GK1
Alternative Catalog Number	ROC-W09-001-GK1
Manufacturer	Rockland Immunochemicals
Category	Zellen/Zellkultur
Conjugation	Unconjugated
Concentration	4.0 mg/mL by modified Lowry assay
Buffer	1X RIPA Buffer with HALT Protease and Phosphatase Inhibitors
Source	Human
Form	Liquid (sterile filtered)
Application Dilute	ChIP: User Optimized, IP: User Optimized, WB: User Optimized
Application Notes	Multi-purpose LNCaP nuclear extracts are especially prepared as positive control for multiple assays including western blot, immunoprecipitation (IP), capture ELISA or other assays requiring native protein sample. For separation by SDS-PAGE and subsequen