



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## MOLT-4 Cell Nuclear Extract ROC-W09-001-GK3

Article Name	MOLT-4 Cell Nuclear Extract
Biozol Catalog Number	ROC-W09-001-GK3
Supplier Catalog Number	W09-001-GK3
Alternative Catalog Number	ROC-W09-001-GK3
Manufacturer	Rockland Immunochemicals
Category	Zellen/Zellkultur
Application	SDS-PAGE, WB
Conjugation	Unconjugated
Concentration	4.0 mg/mL by modified Lowry assay
Buffer	1X RIPA Buffer with HALT Protease and Phosphatase Inhibitors
Source	Human
Form	Liquid (sterile filtered)
Application Dilute	ChIP: User Optimized, IP: User Optimized, WB: User Optimized
Application Notes	W09-001-GK3 has been tested by SDS-PAGE and western blot. Multi- purpose MOLT-4 nuclear extracts are especially prepared as positive control for multiple assays including western blot, immunoprecipitation (IP), capture ELISA or other assays requiring nati