

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Ramos Cell Nuclear Extract

ROC-W09-001-GK5

Article Name	Ramos Cell Nuclear Extract
Biozol Catalog Number	ROC-W09-001-GK5
Supplier Catalog Number	W09-001-GK5
Alternative Catalog Number	ROC-W09-001-GK5
Manufacturer	Rockland Immunochemicals
Category	Zellen/Zellkultur
Conjugation	Unconjugated
Concentration	4.0 mg/mL by modified Lowry assay
Buffer	1X RIPA Buffer with HALT Protease and Phosphatase Inhibitors
Source	Human
Form	Liquid (sterile filtered)
Application Dilute	ChIP: User Optimized, IP: User Optimized, WB: User Optimized
Application Notes	Multi-purpose Ramos nuclear extracts are especially prepared as positive control for multiple assays including western blot, immunoprecipitation (IP), capture ELISA or other assays requiring native protein sample. For separation by SDS-PAGE and subsequen