

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

APOBEC3F, NT (APOBEC3F, DNA dC->dU-editing enzyme APOBEC-3F, Apolipoprotein B mRNA-editing enzyme catalytic polypeptide-like 3F) (MaxLight 650), Rabbit USB-031983-ML650

| Article Name | APOBEC3F, NT (APOBEC3F, DNA dC->dU-editing enzyme APOBEC-3F, Apolipoprotein B mRNA-editing enzyme catalytic polypeptide-like 3F) (MaxLight 650), Rabbit |
|----------------------------|--|
| Biozol Catalog Number | USB-031983-ML650 |
| Supplier Catalog Number | 031983-ML650 |
| Alternative Catalog Number | USB-031983-ML650-100 |
| Manufacturer | US Biological |
| Host | Rabbit |
| Category | Antikörper |
| Application | FLISA, IHC, WB |
| Immunogen | KLH-conjugated synthetic peptide mapping to a fragment of residues within amino acids 22-48 in the N-terminal region of human APOBEC3F. UniProt Accession Q8IUX4 |
| Product Description | MaxLight(TM)650 is a new Far-IR stable dye conjugate comparable to Alexa Fluor(TM)647, DyLight(TM)649, Cy5(TM) and offers better labeling efficiency, brighter imaging and increased immunodetection. Absorbance (655nm), Emission (676nm), Extinction Coe |
| NCBI | 8 |
| UniProt | Q8IUX4 |

| Purity | Purified by Protein A affinity chromatography. |
|--------|---|
| Form | Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight(TM)650. |