

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

A431 cells stimulated with EGF, Unconjugated BYT-ORB198380

Article Name	A431 cells stimulated with EGF, Unconjugated
Biozol Catalog Number	BYT-ORB198380
Supplier Catalog Number	orb198380
Alternative Catalog Number	BYT-ORB198380-100
Manufacturer	Biorbyt
Category	Sonstiges
Application	WB
Conjugation	Unconjugated
Product Description	A431 cells stimulated with EGF...
Concentration	See vial for concentration
Buffer	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Form	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Application Notes	Suggested use: A volume of 10-25ul is suggested for immunoblotting.