

Diagnostica Vertrieb GmbH, Leipziger Straße 4 85386 Eching, Germany Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99 E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

3T3 cells, Unconjugated BYT-ORB198382

Article Name	3T3 cells, Unconjugated
Biozol Catalog Number	BYT-ORB198382
Supplier Catalog Number	orb198382
Alternative Catalog Number	BYT-ORB198382-100
Manufacturer	Biorbyt
Category	Sonstiges
Application	WB
Conjugation	Unconjugated
Product Description	3T3 cells
Concentration	See vial for concentration
Buffer	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Form	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Application Notes	Suggested use: A volume of 10-25ul is suggested for immunobloting.