

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

**E-Mail:** info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## HeLa Cell Nuclear Extract Nocodozole Stimulated, Unconjugated BYT-ORB348686

Article Name	HeLa Cell Nuclear Extract Nocodozole Stimulated, Unconjugated
Biozol Catalog Number	BYT-ORB348686
Supplier Catalog Number	orb348686
Alternative Catalog Number	BYT-ORB348686-200
Manufacturer	Biorbyt
Category	Sonstiges
Application	SDS-PAGE, WB
Conjugation	Unconjugated
Product Description	HeLa Cell Nuclear Extract Nocodozole Stimulated
Concentration	1.0 mg/ml
Buffer	Preservative: None. Stabilizer: 10% (v/v) Glycerol. 1X SDS-PAGE Sample Buffer (62.5 mM Tris HCl, 2% SDS, 10% Glycerol and 0.005% bromophenol blue, pH 6.8)
Source	Human
Purity	The cells were grown in DMEM supplemented with 10% FBS (Fetal Bovine Serum). Cells were treated with 0.2 µg/ml Nocodazole for 30 min. The lysate was prepared by first washing the cells in PBS. Washed cells were then incubated on ice in lysis buffer conta
Form	Liquid (sterile filtered)
Application Dilute	WB: User Optimized

## **Application Notes**

ready-to-use nuclear extracts are especially prepared as positive controls for separation by SDS-PAGE and subsequent western blot analysis. Nuclear extracts are supplied in denaturing buffer without dissociating agents. Heat nuclear extract to 95C for 5