

Diagnostica Vertrieb GmbH, Leipziger Straße 4 85386 Eching, Germany Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99 E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

IDH3G Blocking peptide BYT-ORB1442740

IDH3G Blocking peptide
BYT-ORB1442740
orb1442740
BYT-ORB1442740-500
Biorbyt
Proteine/Peptide
IDH3G Blocking peptide
3421
P51553
Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
The synthetic peptide sequence is selected from aa 50-63 of HUMAN IDH3G
IDH3G
Synthetic peptide