

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## VATC Antibody Blocking peptide BYT-ORB1451780

Article Name	VATC Antibody Blocking peptide
Biozol Catalog Number	BYT-ORB1451780
Supplier Catalog Number	orb1451780
Alternative Catalog Number	BYT-ORB1451780-500
Manufacturer	Biorbyt
Category	Proteine/Peptide
Product Description	VATC Antibody Blocking peptide
NCBI	528
UniProt	P21283
Buffer	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Form	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Sequence	The synthetic peptide sequence used to generate the antibody was selected from the N-term region of human VATC. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Target	ATP6V1C1
Application Notes	Synthetic peptide