



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Mobicol Spin Fixed with screw cap, fixed outlet plug and inserted small 10 $\mu m$ filter MOB-M105010S

Article Name	Mobicol Spin Fixed with screw cap, fixed outlet plug and inserted small 10 $\mu m$ filter
Biozol Catalog Number	MOB-M105010S
Supplier Catalog Number	M105010S
Alternative Catalog Number	MOB-M105010S
Manufacturer	MoBiTec
Category	Chromatographie
Product Description	Empty Mini Columns: Mobicol Spin – Speed-Up your Purification with a Microcentrifuge. The Mobicol Spin columns are developed as fast and convenient routine tool for clean-up of biomolecules by spinning. They can be used without matrix or filled with