



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mobicol with immobilized TPCK-Trypsin MOB-P3702

Article Name	Mobicol with immobilized TPCK-Trypsin
Biozol Catalog Number	MOB-P3702
Supplier Catalog Number	P3702
Alternative Catalog Number	MOB-P3702
Manufacturer	MoBiTec
Category	Chromatographie
Product Description	Enzyme immobilization provides enzyme reutilization and prevents the loss of enzyme activity. The technique further results in increased activity by providing a more suitable microenvironment for the enzyme substrate reaction. Moreover, immobilized s