



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

sequencing primer 2897 (19 mer), lyophilized MOB-PS2897

Article Name	sequencing primer 2897 (19 mer), lyophilized
Biozol Catalog Number	MOB-PS2897
Supplier Catalog Number	PS2897
Alternative Catalog Number	MOB-PS2897
Manufacturer	MoBiTec
Category	Molekularbiologie
Product Description	Phagemid Display System pSKAN - Construction of Libraries with Huge Molecular Diversity. The pSKAN Phagemid Display System allows selection and production of novel proteins with the desired affinities to a given epitope. pSKAN is a combinatorial phag