

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Gel-based E. coli DNA Topoisomerase II (Gyrase) Assay Kit Plus-100 (for) PFD-GDSA100KE

Article Name	Gel-based E. coli DNA Topoisomerase II (Gyrase) Assay Kit Plus-100 (for)
Biozol Catalog Number	PFD-GDSA100KE
Supplier Catalog Number	GDSA100KE
Alternative Catalog Number	PFD-GDSA100KE
Manufacturer	ProFoldin
Category	Kits/Assays
Product Description	This product includes 400 ul of $10 \times Buffer$, $105 \text{ ul of } 250 \text{ ug/ml}$ relaxed DNA, $220 \text{ ul of } 10 \text{ mM}$ ATP, $550 \text{ ul of } 5 \times Gel$ loading buffer and $25 \text{ ul of } 2000 \text{ nM}$ E.coli gyrase