



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

PolySUMO Capture & Release Kit UBP-J4410

Article Name	PolySUMO Capture & Release Kit
Biozol Catalog Number	UBP-J4410
Supplier Catalog Number	J4410
Alternative Catalog Number	UBP-J4410-1
Manufacturer	UBPBio, LLC
Category	Kits/Assays
Species Reactivity	Human
Product Description	The PolySUMO Capture & Release Kit is designed to rapidly enrich polysumoylated proteins from cell or tissue extracts within 6 hours according to a published method (reference 1) with modifications. This approach uses the SIM domains of RNF4 to captu
Buffer	20 mM Tris, 150 mM NaCl, 2 mM BME, 10% Glycerol
Source	E. coli recombinant