

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Cattle betaTG(Beta-Thromboglobulin) ELISA Kit EBT-ELK0743

Article Name	Cattle betaTG(Beta-Thromboglobulin) ELISA Kit
Biozol Catalog Number	EBT-ELK0743
Supplier Catalog Number	ELK0743
Alternative Catalog Number	EBT-ELK0743-96,EBT-ELK0743-48,EBT-ELK0743-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Bovine
Concentration	100 ng/mL
Range	1.57-100 ng/mL
Sensitivity	0.66 ng/mL
Samples	Serum, plasma, tissue homogenates and other biological fluids.

Application Notes	<p>Assay Type: Sandwich. Assay length: 3.5h. Research Area: Cytokine,Hematology,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cattle bTG. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cattle bTG. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cattle bTG, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Cattle bTG in the samples is then determined by comparing the OD of the samples to the standard curve</p>
-------------------	---