

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human Bid(BH3 Interacting Domain Death Agonist) ELISA Kit EBT-ELK1703

Article Name	Human Bid(BH3 Interacting Domain Death Agonist) ELISA Kit
Biozol Catalog Number	EBT-ELK1703
Supplier Catalog Number	ELK1703
Alternative Catalog Number	EBT-ELK1703-96,EBT-ELK1703-48,EBT-ELK1703-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Human
Concentration	100 ng/mL
Range	1.57-100 ng/mL
Sensitivity	0.56 ng/mL
UniProt	<a href="#">P55957</a>
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.

Application Notes	<p>Assay Type: Sandwich. Assay length: 3.5h. Research Area: Signal transduction,Apoptosis,Tumor immunity,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Human Bid. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Human Bid. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Human Bid, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Human Bid in the samples is then determined by comparing the OD of the samples to the standard curve</p>
-------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------