

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Monosan Plus (AP) Polymer anti-Mouse, 6 ml/60 tests MON-APP205**

Article Name	Monosan Plus (AP) Polymer anti-Mouse, 6 ml/60 tests
Biozol Catalog Number	MON-APP205
Supplier Catalog Number	MON-APP205
Alternative Catalog Number	MON-APP205
Manufacturer	Monosan
Category	Sonstiges
Application	IF, IHC-Fr, IHC-P
Product Description	The purpose of the immunohistochemical staining is to make tissue and cell antigens visible. The Plus AP Polymer anti-Mouse is a highly sensitive detection reagent intended for use in immunohistochemistry and immunocytochemistry. The enzyme polymer i...