

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Monosan (AP) Polymer anti-Mouse, 1.000 tests MON-APP206

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name | Monosan (AP) Polymer anti-Mouse, 1.000 tests |
| Biozol Catalog Number | MON-APP206 |
| Supplier Catalog Number | MON-APP206 |
| Alternative Catalog Number | MON-APP206 |
| Manufacturer | Monosan |
| Category | Sonstiges |
| Application | IF, IHC-Fr, IHC-P |
| Product Description | The purpose of the immunohistochemical staining is to make tissue and cell antigens visible. The Plus AP Polymer anti-Mouse is a highly sensitive detection reagent intended for use in immunohistochemistry and immunocytochemistry. The enzyme polymer i... |