



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## MOUNTING MEDIA (WATER BASED) MON-APP901

Article Name	MOUNTING MEDIA (WATER BASED)
Biozol Catalog Number	MON-APP901
Supplier Catalog Number	MON-APP901
Alternative Catalog Number	MON-APP901
Manufacturer	Monosan
Category	Biochemikalien
Application	IHC, IHC-Fr
Product Description	The aqua mounting medium is a widely used for coverslipping from aqueous solutions and is non-fluorescing and has an antifade component to increase the viewing time of the specimen. Use Aqua-Poly/Mount with most fluorescent dyes and stains including