

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

RH GDNF (AF) SHD-100-02AF-10UG

Article Name	RH GDNF (AF)
Biozol Catalog Number	SHD-100-02AF-10UG
Supplier Catalog Number	100-02AF-10UG
Alternative Catalog Number	SHD-100-02AF-10UG
Manufacturer	Fujifilm Irvine Scientific
Category	Sonstiges
Application	Bioassay, SDS-PAGE
Species Reactivity	Human
Product Description	Glial cell-derived neurotrophic factor (GDNF) is a neurotrophic factor that signals through a multicomponent receptor system to activate receptor tyrosine kinase RET signaling. GDNF promotes dopamine uptake, prevents motor neuron apoptosis, and suppo
Molecular Weight	15.2/30.4
UniProt	P39905
Buffer	10 mM Sodium Citrate, 100 mM Sodium Chloride, pH 4.0
Source	E.coli
Purity	95
Form	Lyophilized from a 0.2 μm filtered solution

Sequence	MSPDKQMAVL PRRERNRQAA AANPENSRGK GRRGQRGKNR GCVLTAIHLN VTDLGLGYET KEELIFRYCS GSCDAAETTY DKILKNLSRN RRLVSDKVGQ ACCRPIAFDD DLSFLDDNLV YHILRKHSAK RCGCI
Application Notes	Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.