

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

RH GALECTIN-3 SHD-100-154-500UG

Article Name	RH GALECTIN-3
Biozol Catalog Number	SHD-100-154-500UG
Supplier Catalog Number	100-154-500UG
Alternative Catalog Number	SHD-100-154-500UG
Manufacturer	Fujifilm Irvine Scientific
Category	Sonstiges
Application	Bioassay, SDS-PAGE
Species Reactivity	Human
Product Description	Galectin-3 belongs to the lectin family of carbohydrate binding proteins. Galectin-3 is expressed by a wide range of cell types including activated T cells, tumor cells, macrophages, osteoclasts, fibroblasts, and epithelial cells. Galectin-3 has spec...
Molecular Weight	26.2
UniProt	P17931
Buffer	10 mM Sodium Phosphate, 50 mM Sodium Chloride, pH 7.5
Source	E.coli
Purity	95
Form	Lyophilized from a 0.2 µm filtered solution

Sequence	MADNFSLHDA LSGSGNPNPQ GWPGAWGNQP AGAGGYPGAS YPGAYPGQAP PGAYPGQAPP GAYPGAPGAY PGAPAPGVYP GPPSGPGAYP SSGQPSATGA YPATGPYGAP AGPLIVPYNL PLPGGVVPRM LITILGTVKP NANRIALDFQ RGNDVAFHFN PRFNENNRV IVCNTKLDNN WGREERQSVF PFESGKPFI QVLVEPDHFK VAVNDAHLLQ YNHRVKKLNE I
Application Notes	Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.