



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

RH SHH SHD-100-168-25UG

| Article Name | RH SHH |
|----------------------------|--|
| Biozol Catalog Number | SHD-100-168-25UG |
| Supplier Catalog Number | 100-168-25UG |
| Alternative Catalog Number | SHD-100-168-25UG |
| Manufacturer | Fujifilm Irvine Scientific |
| Category | Sonstiges |
| Application | Bioassay, SDS-PAGE |
| Species Reactivity | Human |
| Product Description | Sonic hedgehog (SHH) is a member of a small group of hedgehog secreted proteins that are essential for development in both vertebrates and invertebrates. There are three mammalian hedgehog homologues, sonic, desert, and indian, that signal via the Pa |
| Molecular Weight | 20.2 |
| UniProt | Q15465 |
| Buffer | 10 mM Sodium Phosphate, pH 7.5 |
| Source | E.coli |
| Purity | 95 |
| Form | Lyophilized from a 0.2 µm filtered solution |

| Sequence | MIIGPGRGFG KRRHPKKLTP LAYKQFIPNV AEKTLGASGR YEGKISRNSE RFKELTPNYN PDIIFKDEEN TGADRLMTQR CKDKLNALAI SVMNQWPGVK LRVTEGWDED GHHSEESLHY EGRALDITTS DRDRSKYGML ARLAVEAGFD WVYYESKAHI HCSVKAENSV AAKSGGCFP |
|-------------------|--|
| Application Notes | Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0. |