

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

RH GRO-G SHD-100-170-500UG

Article Name	RH GRO-G
Biozol Catalog Number	SHD-100-170-500UG
Supplier Catalog Number	100-170-500UG
Alternative Catalog Number	SHD-100-170-500UG
Manufacturer	Fujifilm Irvine Scientific
Category	Sonstiges
Application	Bioassay, SDS-PAGE
Species Reactivity	Human
Product Description	Growth regulated protein gamma (GRO-), also called CXCL3, acts through the chemokine receptor CXCR2 to promote monocyte migration and adhesion. GRO- also induces the migration of cerebellar granule neuron precursor cells
Molecular Weight	7.9
UniProt	P19876
Buffer	10 mM Sodium Phosphate, pH 7.5
Source	E.coli
Purity	95
Form	Lyophilized from a 0.2 µm filtered solution
Sequence	ASVVTELRCQ CLQTLQGIHL KNIQSVNVRS PGPHCAQTEV IATLKNGKKA CLNPASPMVQ KIIEKILNKG STN

Application Notes

Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.