

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

RH EOTAXIN-3 SHD-100-179-100UG

Article Name	RH EOTAXIN-3
Biozol Catalog Number	SHD-100-179-100UG
Supplier Catalog Number	100-179-100UG
Alternative Catalog Number	SHD-100-179-100UG
Manufacturer	Fujifilm Irvine Scientific
Category	Sonstiges
Application	Bioassay, SDS-PAGE
Species Reactivity	Human
Product Description	Eotaxin-3, also known as CCL26, MIP-4-alpha, and TSC-1, is a chemokine that is made by vascular endothelial and lung epithelial cells following interleukin 4 (IL-4) or interleukin 13 (IL-13) stimulation. Eotaxin-3 signals through the G protein-couple
Molecular Weight	8.4
UniProt	Q9Y258
Buffer	0.1% Trifluoroacetic Acid (TFA)
Source	E.coli
Purity	95
Form	Lyophilized from a 0.2 µm filtered solution
Sequence	TRGSDISKTC CFQYSHKPLP WTWVRSYEFT SNSCSQRAVI FTTKRGKKVC THPRKKWVQK YISLLKTPKQ L

Application Notes

Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.