

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

RH BAFF (AF) SHD-100-24AF-1MG

| | |
|----------------------------|---|
| Article Name | RH BAFF (AF) |
| Biozol Catalog Number | SHD-100-24AF-1MG |
| Supplier Catalog Number | 100-24AF-1MG |
| Alternative Catalog Number | SHD-100-24AF-1MG |
| Manufacturer | Fujifilm Irvine Scientific |
| Category | Sonstiges |
| Application | Bioassay, SDS-PAGE |
| Species Reactivity | Human |
| Product Description | B cell-activating factor (BAFF), or B lymphocyte stimulator (BLyS), is a type II member of the tumor necrosis factor (TNF) superfamily. BAFF is expressed as a transmembrane protein on T cells, macrophages, and dendritic cells. The transmembrane domai... |
| Molecular Weight | 18.5 |
| UniProt | Q9Y275 |
| Buffer | 10 mM Sodium Phosphate, pH 7.5 |
| Source | E.coli |
| Purity | 90 |
| Form | Lyophilized from a 0.2 µm filtered solution |

| | |
|-------------------|--|
| Sequence | MHHHHHHLVP RAVQGPEETV TQDCLQLIAD SETPTIQKGS YTFVPWLLSF KRGSALEEKE NKILVKETGY FFIYGQVLYT DKTYAMGHLI QRKKVHVFGD ELSLVTLFRC IQNMPETLPN NSCYSAGIAK LEEGDELQLA IPRENAQISL DGDVTFFGAL KLL |
| Application Notes | Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0. |