

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

RH BAFF (AF) SHD-100-24AF-500UG

Article Name	RH BAFF (AF)
Biozol Catalog Number	SHD-100-24AF-500UG
Supplier Catalog Number	100-24AF-500UG
Alternative Catalog Number	SHD-100-24AF-500UG
Manufacturer	Fujifilm Irvine Scientific
Category	Sonstiges
Application	Bioassay, SDS-PAGE
Species Reactivity	Human
Product Description	B cell-activating factor (BAFF), or B lymphocyte stimulator (BLyS), is a type II member of the tumor necrosis factor (TNF) superfamily. BAFF is expressed as a transmembrane protein on T cells, macrophages, and dendritic cells. The transmembrane domai...
Molecular Weight	18.5
UniProt	Q9Y275
Buffer	10 mM Sodium Phosphate, pH 7.5
Source	E.coli
Purity	90
Form	Lyophilized from a 0.2 µm filtered solution

Sequence	MHHHHHHLVP RAVQGPEETV TQDCLQLIAD SETPTIQKGS YTFVPWLLSF KRGSALEEKE NKILVKETGY FFIYGQVLYT DKTYAMGHLI QRKKVHVFGD ELSLVTLFRC IQNMPETLPN NSCYSAGIAK LEEGDELQLA IPRENAQISL DGDVTFFGAL KLL
Application Notes	Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.