

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

RH GALECTIN-1 (AF) SHD-100-30AF-100UG

Article Name	RH GALECTIN-1 (AF)
Biozol Catalog Number	SHD-100-30AF-100UG
Supplier Catalog Number	100-30AF-100UG
Alternative Catalog Number	SHD-100-30AF-100UG
Manufacturer	Fujifilm Irvine Scientific
Category	Sonstiges
Application	Bioassay, SDS-PAGE
Species Reactivity	Human
Product Description	Galectin-1 belongs to the lectin family of carbohydrate binding proteins and binds glycans as both a monomer and a homodimer. Galectin-1 is produced in peripheral lymphoid organs and inflammatory sites. Galectin-1 plays important roles in acute and c...
Molecular Weight	14.7
UniProt	P09382
Buffer	10 mM Sodium Phosphate, pH 7.5
Source	E.coli
Purity	95
Form	Lyophilized from a 0.2 µm filtered solution

Sequence	MACGLVASNL NLKPGECLRV RGEVAPDAKS FVLNLGKDSN NLCLHFNPRF NAHGDANTIV CNSKDGGAWG TEQREAVFPF QPGSVAEVC TFDQANLTVK LPDGYEFKFP NRLNLEAINY MAADGDFKIK CVAFD
Application Notes	Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.