

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### RM RELM-A SHD-200-28-1MG

|                            |   |
|----------------------------|---|
| Article Name               | RM RELM-A   |
| Biozol Catalog Number      | SHD-200-28-1MG  |
| Supplier Catalog Number    | 200-28-1MG  |
| Alternative Catalog Number | SHD-200-28-1MG  |
| Manufacturer               | Fujifilm Irvine Scientific  |
| Category                   | Sonstiges   |
| Application                | Bioassay, SDS-PAGE  |
| Species Reactivity         | Mouse   |
| Product Description        | Resistin-like molecule-alpha (RELM-alpha) is a member of the RELM family of secreted proteins containing conserved C-terminus cysteines. The RELM family consists of Resistin (FIZZ3), RELM-alpha (FIZZ1), RELM-beta (FIZZ2), and RELM-gamma (FIZZ4). RELM... |
| Molecular Weight           | 9.6   |
| UniProt                    | <a href="#">Q9EP95</a>  |
| Buffer                     | 10 mM Sodium Phosphate, pH 7.5  |
| Source                     | E.coli  |
| Purity                     | 95  |
| Form                       | Lyophilized from a 0.2 µm filtered solution   |
| Sequence                   | MDETIEIIVE NKVKELLANP ANYPSTVTKT LSCTSVKTMN RWASCPAGMT<br>ATGCACGFAC GSWEIQSGDT CNCLCLLV DW T TARCCQLS  |

|                   |  |
|-------------------|--|
| Application Notes | Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0. |
|-------------------|--|