

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

RM IL-36 G (AF) SHD-200-81AF-100UG

| Article Name | RM IL-36 G (AF) |
|----------------------------|--|
| Biozol Catalog Number | SHD-200-81AF-100UG |
| Supplier Catalog Number | 200-81AF-100UG |
| Alternative Catalog Number | SHD-200-81AF-100UG |
| Manufacturer | Fujifilm Irvine Scientific |
| Category | Proteine/Peptide |
| Application | Bioassay, SDS-PAGE |
| Species Reactivity | Mouse |
| Product Description | Interleukin 36 gamma (IL-36) is a member of the interleukin 1 (IL-1) cytokine family and protects against pathogens in the skin, lung, and stomach epithelial barriers. IL-36 binds the interleukin-1 receptor accessory protein (IL-1RAcP) and the orph |
| Molecular Weight | 17.5 |
| UniProt | Q8R460 |
| Buffer | 0.1% Trifluoroacetic Acid (TFA) |
| Source | E.coli |
| Purity | 95 |
| Form | Lyophilized from a 0.2 µm filtered solution |

| Sequence | MGRETPDFGE VFDLDQQVWI FRNQALVTVP RSHRVTPVSV TILPCKYPES LEQDKGIAIY LGIQNPDKCL FCKEVNGHPT LLLKEEKILD LYHHPEPMKP FLFYHTRTGG TSTFESVAFP GHYIASSKTG NPIFLTSKKG EYYNINFNLD IKS |
|-------------------|--|
| Application Notes | Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0. |