

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### B2M Peptide PRS-43-039P

Article Name	B2M Peptide
Biozol Catalog Number	PRS-43-039P
Supplier Catalog Number	43-039P
Alternative Catalog Number	PRS-43-039P-0.1
Manufacturer	ProSci
Category	Proteine/Peptide
Application	Blocking
NCBI	<a href="#">567</a>
UniProt	<a href="#">P61769</a>
Buffer	B2M Peptide is supplied as 100 µg of dried peptide. When the peptide is reconstituted in 200 µL water, the concentration will be 0.5 mg/mL.
Form	Lyophilized
Application Notes	B2M Antibody ships at 0.5 mg/mL and the B2M Peptide ships as a 100 µg pellet. When peptide is reconstituted in 200 µL water, the concentration will be 0.5 mg/mL, identical to that of the antibody. The best starting ratio would be 1:1 (which means molar excess of peptide relative to antibody when identical volumes are mixed). Mix equal volume of peptide and antibody at the required dilution and leave at ambient temperature for 20 minutes. Then incubate and develop in parallel two blots: one with the antibody and one with the antibody which has been pre-absorbed with the provided peptide.