



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## GAPDH Peptide PRS-45-061P

Article Name	GAPDH Peptide
Biozol Catalog Number	PRS-45-061P
Supplier Catalog Number	45-061P
Alternative Catalog Number	PRS-45-061P-0.1
Manufacturer	ProSci
Category	Proteine/Peptide
Application	Blocking
NCBI	2597
UniProt	P04406
Buffer	GAPDH peptide is supplied as 100 $\mu g$ of dried peptide. When the peptide is reconstituted in 200 $\mu L$ water, the concentration will be 0.5 mg/mL.
Form	Lyophilized

**Application Notes** 

GAPDH Peptide is used for blocking the activity of GAPDH antibody (catalog number 45-061). The GAPDH antibody ships at 0.5 mg/mL and the GAPDH Peptide ships as a 100 ug pellet. When peptide is reconstituted in 200 uL water, the concentration will be 0.5 mg/mL, identical to that of the antibody. The best starting ratio would be 1:1 (which means molar excess of peptide relative to antibody when identical volumes are mixed). Mix equal volume of peptide and antibody at the required dilution and leave at ambient temperature for 20 minutes. Then incubate and develop in parallel two blots: one with the antibody and one with the antibody which has been preabsorbed with the provided peptide.