



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## HAP1 Peptide PRS-45-722P

Article Name	HAP1 Peptide
Biozol Catalog Number	PRS-45-722P
Supplier Catalog Number	45-722P
Alternative Catalog Number	PRS-45-722P-0.1
Manufacturer	ProSci
Category	Proteine/Peptide
Application	Blocking
NCBI	15114
UniProt	O35668
Buffer	HAP1 peptide is supplied as 100 $\mu g$ of dried peptide. When the peptide is reconstituted in 200 $\mu L$ water, the concentration will be 0.5 mg/mL.
Form	Lyophilized

**Application Notes** 

HAP1 Peptide is used for blocking the activity of HAP1 antibody (catalog number 45-722). The HAP1 antibody ships at 0.5 mg/mL and the HAP1 Peptide ships as a 100 ug pellet. When peptide is reconstituted in 200 uL water, the concentration will be 0.5 mg/mL, identical to that of the antibody. The best starting ratio would be 1:1 (which means molar excess of peptide relative to antibody when identical volumes are mixed). Mix equal volume of peptide and antibody at the required dilution and leave at ambient temperature for 20 minutes. Then incubate and develop in parallel two blots: one with the antibody and one with the antibody which has been preabsorbed with the provided peptide.