

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Phospho Serine Antibody, Unconjugated, Rabbit, Polyclonal PRS-54-207

Article Name	Phospho Serine Antibody, Unconjugated, Rabbit, Polyclonal
Biozol Catalog Number	PRS-54-207
Supplier Catalog Number	54-207
Alternative Catalog Number	PRS-54-207-0.1
Manufacturer	ProSci
Host	Rabbit
Category	Antikörper
Application	ELISA, ICC, IHC, WB
Species Reactivity	Human, Mouse
Immunogen	Phosphoserine-KLH and phosvitin mixture
Conjugation	Unconjugated
Clonality	Polyclonal
Concentration	batch dependent
Buffer	PBS Buffer containing 0.02% sodium azide.
Form	Liquid
Application Dilute	Optimal dilutions for each application to be determined by the researcher.

Application Notes

Phospho-Serine antibody can be used at 1:500 to 1:1,000 for ELISA, 1:200-1:1,000 for Western Blot testing and 1ug/mg sample for IP. Customer using ECL system may need to further dilute the Abs. User needs to try out the titration for optimal concentration.