

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### ULK1 Peptide PRS-7335P

|                            |                                                                                                                                                                                                                                       |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | ULK1 Peptide                                                                                                                                                                                                                          |
| Biozol Catalog Number      | PRS-7335P                                                                                                                                                                                                                             |
| Supplier Catalog Number    | 7335P                                                                                                                                                                                                                                 |
| Alternative Catalog Number | PRS-7335P-0.05                                                                                                                                                                                                                        |
| Manufacturer               | ProSci                                                                                                                                                                                                                                |
| Category                   | Proteine/Peptide                                                                                                                                                                                                                      |
| Application                | Blocking                                                                                                                                                                                                                              |
| Concentration              | 200 ug/mL                                                                                                                                                                                                                             |
| NCBI                       | <a href="#">8408</a>                                                                                                                                                                                                                  |
| UniProt                    | <a href="#">O75385</a>                                                                                                                                                                                                                |
| Buffer                     | PBS pH 7.2 (10 mM NaH <sub>2</sub> PO <sub>4</sub> , 10 mM Na <sub>2</sub> HPO <sub>4</sub> , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide                                                                |
| Form                       | Liquid                                                                                                                                                                                                                                |
| Application Notes          | ULK1 peptide is used for blocking the activity of ULK1 antibody (catalog number 7335). It usually blocks the antibody activity completely in Western blot by incubating the peptide with equal volume of antibody for 30 min at 37°C. |