

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## **Eotaxin-3 Antibody (biotin), Biotin, Goat, Polyclonal PRS-XP-5127BT**

| Article Name               | Eotaxin-3 Antibody (biotin), Biotin, Goat, Polyclonal   |
|----------------------------|---|
| Biozol Catalog Number      | PRS-XP-5127BT   |
| Supplier Catalog Number    | XP-5127Bt   |
| Alternative Catalog Number | PRS-XP-5127BT-0.05  |
| Manufacturer               | ProSci  |
| Host                       | Goat  |
| Category                   | Antikörper  |
| Application                | ELISA, WB   |
| Species Reactivity         | Human   |
| Immunogen                  | Produced from sera of goats pre-immunized with highly pure (>98%) recombinant hEotaxin-3 (human Eotaxin-3). |
| Conjugation                | Biotin  |
| Clonality                  | Polyclonal  |
| Concentration              | batch dependent   |
| NCBI                       | 10344   |
| UniProt                    | Q9Y258  |
| Form                       | Lyophilized   |
| Application Dilute         | Centrifuge vial prior to opening.   |

**Application Notes** 

ELISA:Sandwich: To detect hEotaxin-3 by sandwich ELISA (using 100  $\mu\text{L/well}$  antibody solution) a concentration of 0.25 - 1.0  $\mu\text{g/mL}$  of this antibody is required. This biotinylated polyclonal antibody, in conjunction with our Polyclonal Anti-Human Eotoax