

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Fas Ligand Antibody (biotin), Biotin, Goat, Polyclonal PRS-XP-5130BT

|                            |   |
|----------------------------|---|
| Article Name               | Fas Ligand Antibody (biotin), Biotin, Goat, Polyclonal  |
| Biozol Catalog Number      | PRS-XP-5130BT   |
| Supplier Catalog Number    | XP-5130Bt   |
| Alternative Catalog Number | PRS-XP-5130BT-0.05  |
| Manufacturer               | ProSci  |
| Host                       | Goat  |
| Category                   | Antikörper  |
| Application                | ELISA, WB   |
| Species Reactivity         | Human   |
| Immunogen                  | Produced from sera of goats pre-immunized with highly pure (>98%) recombinant hFasL/Apo1 (human FasL/Apo1). |
| Conjugation                | Biotin  |
| Clonality                  | Polyclonal  |
| Concentration              | batch dependent   |
| NCBI                       | <a href="#">356</a>   |
| UniProt                    | <a href="#">P48023</a>  |
| Form                       | Lyophilized   |
| Application Dilute         | Centrifuge vial prior to opening.   |

|                   |   |
|-------------------|---|
| Application Notes | <p>ELISA:Sandwich: To detect hFasL/Apo1L by sandwich ELISA (using 100 µL/well antibody solution) a concentration of 0.25 - 1.0 µg/mL of this antibody is required. This biotinylated polyclonal antibody, in conjunction with our Polyclonal Anti-Human FasL/Apo1L (XP-5130) as a capture antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hFasL/Apo1L.</p> <p>Western Blot:To detect hFasL/Apo1 by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 µg/mL. Used in conjunction with compatible secondary reagents the detection limit for recombinant hFasL/Apo1 is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions.</p> |
|-------------------|---|