

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Control Normal Ilama IgG, Whole Molecule (FITC), Camelus PRS-XR-8005

Article Name	Control Normal llama IgG, Whole Molecule (FITC), Camelus
Biozol Catalog Number	PRS-XR-8005
Supplier Catalog Number	XR-8005
Alternative Catalog Number	PRS-XR-8005-0.1
Manufacturer	ProSci
Host	Camelus
Category	Antikörper
Application	ELISA, FC, IF, WB
Concentration	1 mg/mL
Molecular Weight	547.7
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 0.01% (w/v) Sodium Azide
Form	Lyophilized
Application Notes	Llama IgG Fluorescein can be utilized as a control or standard reagent in SDS, Western Blotting, Flow Cytometry, Immunofluorescence, and FLISA experiments. Specific conditions for reactivity should be optimized by the end user.