



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

ClickTech Oligo Link Kit L BCL-BCK-OL-L

Article Name	ClickTech Oligo Link Kit L
Biozol Catalog Number	BCL-BCK-OL-L
Supplier Catalog Number	BCK-OL-L
Alternative Catalog Number	BCL-BCK-OL-L
Manufacturer	baseclick
Category	Kits/Assays
Product Description	This ClickTech Oligo Link Kit can be used for e.g.: a) click-labeling of Oligonucleotides. b) ligating two Oligonucleotides. Compared to the enzymatic ligation "click ligation exhibits several advantages: no byproducts, increased ligated product and
Molecular Weight	n.a.
Form	n.a.
CAS Number	[n.a.]
Formula	n.a.
Application Notes	Modification of Oligonucleotides