



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Anti-Kaede pAb MBL-PM012M

Article Name	Anti-Kaede pAb
Biozol Catalog Number	MBL-PM012M
Supplier Catalog Number	PM012M
Alternative Catalog Number	MBL-PM012M
Manufacturer	MBL
Host	Rb
Category	Antikörper
Application	WB
Immunogen	Full-length Kaede
Product Description	Anti-Kaede pAb, Rabbit, Polyclonal
Clonality	Polyclonal
Isotype	IgG
Buffer	$100~\mu L$ volume of $10~mM$ phosphate buffer containing 0.3 M NaCl and 0.09 % NaN3, pH 8.0. *Azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush plenty of water when disposing materials containing azi
Source	This antibody was fractioned by salting out and gel filtration on a Sephadex G-200 column. The rabbit was immunized with the recombinant CoralHue Kaede.
Target	Kaede

Application Dilute

WB: 1:1000 (chemiluminesence detection system)