



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

Anti-DDDDK-tag pAb MBL-PM020

| Article Name               | Anti-DDDDK-tag pAb   |
|----------------------------|--|
| Biozol Catalog Number      | MBL-PM020  |
| Supplier Catalog Number    | PM020  |
| Alternative Catalog Number | MBL-PM020  |
| Manufacturer               | MBL  |
| Host                       | Rb   |
| Category                   | Antikörper   |
| Application                | ICC, IP, WB  |
| Immunogen                  | Synthetic peptide binding KLH, DYKDDDDK  |
| Product Description        | Anti-DDDDK-tag pAb, Rabbit, Polyclonal   |
| Clonality                  | Polyclonal   |
| Isotype                    | IgG  |
| Buffer                     | $100\;\mu\text{l}$ volume of PBS containing 50% glycerol, pH 7.2. Contains no preservative.  |
| Source                     | This antibody was purified from rabbit serum using affinity column. The rabbit was immunized with carrier protein (CP) conjugated synthetic peptide, XXXDDDDK-CP, CP-XXXDDDDK. |
| Target                     | DDDDK  |
| Application Dilute         | WB: 1:1,000 (chemiluminescence detection system). IP: 5 $\mu\text{L/sample}.$ ICC: 1:1000  |