

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### RNA Marker Low ABN-R0001

Article Name	RNA Marker Low
Biozol Catalog Number	ABN-R0001
Supplier Catalog Number	R0001
Alternative Catalog Number	ABN-R0001-50
Manufacturer	Abnova
Category	Molekularbiologie
Application	EP
Product Description	The RNA Marker Low consists of seven single-stranded RNAs. The 20-base and 50-base RNA are synthesized by chemically (not phosphorylated). The 100, 200, 300, 400 and 500 bases are synthesized by in vitro transcription. The RNA Marker Low is suitable ...
Buffer	In 10 mM Tris-HCl (pH 8.0) buffer containing 1 mM EDTA.
Application Dilute	1 uL (0.7 ug) RNA marker low mix with 5 uL of gel loading buffer.