



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Protein Staining Solution ABN-TM0008

Article Name	Protein Staining Solution
Biozol Catalog Number	ABN-TM0008
Supplier Catalog Number	TM0008
Alternative Catalog Number	ABN-TM0008-500
Manufacturer	Abnova
Category	Sonstiges
Application	EP
Product Description	Protein Staining Solution stains proteins in polyacrylamide gel after SDS gel electrophoresis (detection limit >= 8 ng of protein) with *CoomassieR-G250. All processes including washing and destaining processes can be performed in 1.5 hr. After stain